

# Manual Procedure

<b>Cat. No. 12741</b>	R1	1 x	10 ml
<b>For 90 tests</b>	R2	6 x Powder for	1.5 ml
	R3	1 x	7 ml

## Copper in Urine

Without deproteinization

3,5-diBr.-PAESA Colorimetric method

### Powder Reagents

#### Test principle

At pH 4.70, copper is released from the ceruloplasmin complex. Ascorbic acid acts to reduce the cupric copper to cuprous state which reacts with 3,5-diBr.-PAESA to produce a colored complex.

#### Concentrations in the test

<b>Reagent R1</b>			
Acetate buffer pH = 5.7	60	mmol/L	
Detergent			
<b>Reagent R2</b>			
Ascorbic acid	50	mmol/L	
<b>Reagent R3</b>			
4-(3,5-diBrom-2-pyridyl)-N-sulfopropylanilin (PAESA)	0.18	mmol/L	
<b>Standard:</b> Concentration: As indicated on the bottle (5 µg/dl)			

#### Stability and preparation of working reagent

**Reagent R1:** liquid.

**Reagent R2:** powder.

**Reagent R3:** liquid ready to use.

All reagents are stable up to expiry date given on label when stored at +20 → +25 °C.

#### Working Reagent:

Add 1.5 ml from bottle R1(Buffer) into one vial R2 (Powder). and mix. Stability: 7 days at 20 - 25 °C.

#### Specimen collection and handling

Urine: 24 hours urine collected in plastic container, preferably polyethylene.

Centrifuge and use the supernatant for the test.

Stability: 24 hour at 2 – 8 °C. 15 days at -20°C.

#### Standard

Copper STD. Cat. No. 16081

#### Procedure

Wavelength	Hg 578 (570 - 590 nm)
Spectrophotometer	582nm
Cuvette	1 cm light path
Temperature	20 – 25 °C
Measurement	against reagent blank
Reaction	End point

#### Assay

	Blank	Standard	Sample Blank	Sample
Distilled water	1000 µl	--	100 µl	--
Standard	--	1000 µl	--	--
Sample	--	--	1000 µl	1000 µl
Reagent R3	100 µl	100 µl	--	100 µl
Working Reagent	100 µl	100 µl	100 µl	100 µl

Mix, incubate for 5 min. at 20 – 25 °C. Read the absorbance of sample blank (A<sub>SB</sub>) against distilled water, and the absorbance of sample (A<sub>S</sub>) and standard (A<sub>ST</sub>) against reagent Blank. The final colour is stable for at least 30 min.

#### Calculation

$$\text{Conc. Copper, Urine (µg/dl)} = \frac{A_S - A_{SB}}{A_{ST}} \times \text{Conc. STD. /Urine (µg/dl)}$$

$$\text{Conc. Copper, Urine/ 24 hr.} = \frac{\text{Copper µg/dl (Urine)} \times (\text{vol./ ml})_{\text{urine/ 24 hr.}}}{100} \text{ µg/24 hr.}$$

#### Linearity

Up to 20 µg/dl.

If the result exceeds 20 µg/dl, repeat the test using diluted sample (1+1) with sodium chloride solution (0.9 %) and multiply the result by 2.

#### Precautions

Please use copper free reaction tubes and cuvettes.

#### Interference

- Bilirubin up to 5 mg/dl has no effect on this procedure. Use sample blank for higher bilirubin concentration.
- Lipemic serum should have a sample blank for each test.

#### Reference range

Urine/ 24 hr.	10 – 60	µg/24 hr.
---------------	---------	-----------

#### References

- Abe, A., Yamashita S., and al., Sensitive, Direct Colorimetric Assay for Copper in Serum, Clin. Chem., 35, (1989). 552.
- Landers JW, Zak B. Determination of serum copper and iron in a single small sample. Amer J Clin Path 1958, 29: 590 - 2.
- Houwen RHJ, Hattun van I, Hoogenraad TU. Wilson disease. Netherlands J. Med 1993, 43: 26 - 37.
- Tanzi RE, et al. The Wilson disease gene is a copper transporting ATPase with homology to the Menkes disease gene. Nature Genetics 1993, 5: 344 - 50.
- Aggett PT. Aspects of neonatal metabolism of trace elements. Acta Paediatr.